


**VETERINARY DENTISTRY**  
**Cats & Dogs**  
*Catalogue*





*Take care opening the patient's mouth - too wide may cause damage*

**DRANBS Mouth Gag**

Art. No	Description
282410	DRANBS mouth gag, slim model, 5 cm
282411	DRANBS mouth gag, slim model, 8 cm
282412	DRANBS mouth gag, slim model, 12 cm
282413	DRANBS mouth gag, strong model, 7 cm
282414	DRANBS mouth gag, strong model, 9 cm
282415	DRANBS mouth gag, strong model, 12 cm
282416	DRANBS mouth gag, strong model, 17 cm




**DRANBS Lip Retractor with Soft Coated Jaws**

Art. No	Description
282417	DRANBS lip retractor with soft coated jaws




**DRANBS Minnesota Retractor**

Art. No	Description
282385	DRANBS Minnesota retractor




DRANBS Dental Mirror

Art. No	Description
282402	DRANBS dental mirror, model Fahrenkrug, dual ended



DRANBS Explorer and Measuring Probe

Art. No	Description
<b>C</b> 282407	DRANBS explorer and measuring probe, extra long, 17mm
<b>C</b> 282380	DRANBS explorer and measuring probe, 10mm
<b>F</b> 282395	DRANBS explorer and measuring probe, fine 10mm



DRANBS Resorptive Lesion Probe

Art. No	Description
282397	DRANBS resorptive lesion probe



DRANBS Sickle Scaler, Curved

Art. No	Description
<b>C</b> 282382	DRANBS Sickle scaler, curved
<b>F</b> 282398	DRANBS Sickle scaler, curved, fine, cats/small dogs



DRANBS Jaquette Scaler

Art. No	Description
282383	DRANBS Jaquette scaler, straight blade, supragingival



DRANBS Kaplan Scalers

Art. No	Description
282401	DRANBS Kaplan scaler, dual ended, for large patients like dogs



DRANBS Tooth Scaler

Art. No	Description
282403	DRANBS tooth scaler, right
282404	DRANBS tooth scaler, left



DRANBS McCall Scaler

Art. No	Description
282406	DRANBS McCall scaler, dual ended



DRANBS Curette

Art. No	Description
282384	DRANBS standard Jaquette curette 31/32, subgingival
282399	DRANBS Gracey curette 13/14, subgingival, fine, cats/small dogs



DRANBS Tartar Removing Forceps

Art. No	Description
282386	DRANBS tartar removing forceps



DRANBS Periosteal Elevator Double Ended

Art. No	Description
282381	DRANBS periosteal elevator double ended
282396	DRANBS periosteal elevator, double ended, fine, cats/small dogs



DRANBS Periosteal Elevator, Molt, No 9

Art. No	Description
282390	DRANBS periosteal elevator, Molt, No 9



DRANBS Periosteal Elevator, Rougine, double ended, Nordberg



Art. No	Description
282405	DRANBS periosteal elevator, Rougine, double ended, Nordberg

DRANBS Periosteal Elevator, Standard Handle

Art. No	Description
282355	DRANBS periosteal elevator, standard handle



DRANBS Deciduous Elevators, Curved




Art. No	Description
282391	DRANBS curved deciduous elevator, S
282392	DRANBS curved deciduous elevator, M
282393	DRANBS curved deciduous elevator, L



### Single Ended Round Elevators

Every vet working with dentistry has a favourite shape of the handles of elevators and luxators. The most important considerations are a good ergonomic fit to secure a firm grip and the sharpness of the tip. When choosing the DRANBS pear shaped handle, you get a perfect fit in the palm of your hand as well as a choice between two handle lengths. This will allow you to gain complete control of the tip of the instrument and you will get a good experience of holding the optimal tooth extraction instrument in your hand.



### DRANBS Elevator, Winged, Standard Handle

Art. No	Description
282330	DRANBS elevator, winged, 1 mm, standard handle
282331	DRANBS elevator, winged, 2 mm, standard handle
282332	DRANBS elevator, winged, 3 mm, standard handle
282333	DRANBS elevator, winged, 4 mm, standard handle
282334	DRANBS elevator, winged, 5 mm, standard handle
282335	DRANBS elevator, winged, 6 mm, standard handle

### Sharpening Stone



With a flat surface, this dental sharpening stone can be used for sharpening most dental hand instruments including curettes, elevators and scalars.

In order for the instruments to maintain optimal function, they should be kept sharp at all times. For best results, sharpen the instruments after use. The stone can be used in dry condition or together with water or oil. See page 68 for more sharpening stones.

Art. No	Description
282419	DRANBS sharpening stone, 15 x 5 x 0.75 cm




DRANBS Elevator, Winged, Stubby Handle

Art. No	Description
282336	DRANBS elevator, winged, 1 mm, stubby handle
282337	DRANBS elevator, winged, 2 mm, stubby handle
282338	DRANBS elevator, winged, 3 mm, stubby handle
282339	DRANBS elevator, winged, 4 mm, stubby handle
282340	DRANBS elevator, winged, 5 mm, stubby handle
282341	DRANBS elevator, winged, 6 mm, stubby handle



DRANBS Elevator Set, Winged, 1- 6 mm, Stubby Handle, in Stainless Steel Case

Art. No	Description
282356	DRANBS elevator set, winged, 1 - 6 mm, stubby handle, in stainless steel case




DRANBS Elevator, Straight Tip, Stubby Handle

Art. No	Description
282347	DRANBS elevator, straight tip, 1 mm, stubby handle
282348	DRANBS elevator, straight tip, 2 mm, stubby handle




DRANBS luxator, stubby handle

Art. No	Description
282342	DRANBS luxator, 1 mm, stubby handle
282343	DRANBS luxator, 2 mm, stubby handle
282344	DRANBS luxator, 3 mm, stubby handle
282345	DRANBS luxator, 4 mm, stubby handle
282346	DRANBS luxator, 5 mm, stubby handle



DRANBS Serrated Luxators, Stubby Handle

Art. No	Description
282324	DRANBS serrated luxator, 1 mm, stubby handle
282325	DRANBS serrated luxator, 2 mm, stubby handle
282326	DRANBS serrated luxator, 3 mm, stubby handle
282327	DRANBS serrated luxator, 4 mm, stubby handle
282328	DRANBS serrated luxator, 5 mm, stubby handle
282329	DRANBS serrated luxator, 6 mm, stubby handle



DRANBS Straight Root Tip Pick, Standard Handle

Art. No	Description
282354	DRANBS straight root tip pick, standard handle




DRANBS Extraction Forceps

Art. No	Description
<b>C</b> 282387	DRANBS extraction forceps, 13 cm
<b>F</b> 282400	DRANBS extraction forceps, fine, cats/small dogs, 12 cm




DRANBS Extraction Forceps, Right Angled

Art. No	Description
282394	DRANBS extraction forceps, right angled




DRANBS Root Tip Extraction Forceps

Art. No	Description
282388	DRANBS root tip extraction forceps




DRANBS Fragment Extraction Forceps

Art. No	Description
282389	DRANBS fragment extraction forceps




### Rabbits and Rodents

Working with rabbits and rodents demands special hand instruments. DRANBS has developed a range, which covers many of the procedures in this growing field of veterinary dentistry.



### DRANBS Incisor Luxator

Art. No	Description
282300	DRANBS incisor luxator



### DRANBS Molar Luxator

Art. No	Description
282301	DRANBS molar luxator



### DRANBS Spatula

Art. No	Description
282302	DRANBS spatula



DRANBS Cheek Dilators

Art. No	Description
282307	DRANBS cheek dilator, S
282308	DRANBS cheek dilator, M
282309	DRANBS cheek dilator, L



DRANBS Cheek Dilator, Curved

Art. No	Description
282312	DRANBS cheek dilator, curved, medium, 2.2 cm
282310	DRANBS cheek dilator, curved, long, 4.1 cm
282311	DRANBS cheek dilator, curved, x-long, 6.0 cm



DRANBS Mouth Gag

Art. No	Description
282314	DRANBS rodent mouth gag, wire



DRANBS Rodent Mouth Gag

Art. No	Description
282315	DRANBS rodent mouth gag



DRANBS Table Top Mouth Gag

Art. No	Description
282318	DRANBS table top mouth gag, base 2.5 cm



DRANBS Table Top Mouth Gag, small

Art. No	Description
282321	DRANBS table top mouth gag, small, base 1.5 cm



DRANBS Rabbit Molar Cutter

Art. No	Description
282320	DRANBS rabbit molar cutter, 14 cm
282316	DRANBS rabbit molar cutter, 17 cm, L



DRANBS Rabbit Molar Teeth Extraction Forceps

Art. No	Description
282317	DRANBS rabbit molar teeth extraction forceps




282303

DRANBS Rabbit and Rodent Tooth Rasps

Art. No	Description
282303	DRANBS rodent tooth rasp, narrow
282304	DRANBS rodent tooth rasp, wide



282304



### Instrument Sets

Different dental procedures require various assortments of instruments. In DRANBS we have the most common instruments in sets.

The sets are sold in nice boxes. It is advisable to sterilize and store the instruments in a clean environment before and in between use.

Each instrument can be purchased individually.



*Each instrument can be purchased individually*

#### DRANBS Rabbit and Rodent Set, 8 Pieces

The set consists of:

- 282300 DRANBS Incisor Luxator
- 282301 DRANBS Molar Luxator
- 282302 DRANBS Spatula
- 282307 DRANBS Cheek Dilator, S
- 282308 DRANBS Cheek Dilator, M
- 282314 DRANBS rodent mouth gag, wire
- 282315 DRANBS Rodent Mouth Gag
- 282317 DRANBS Rabbit Molar Teeth Extraction Forceps.

Art. No	Description
282319	DRANBS rabbit and rodent set, 8 pieces



*Each instrument can be purchased individually*

#### DRANBS Instrument Set, 14 Pieces, Stubby Handle

The set consist of:

- 282336 - 282340 DRANBS Elevator, winged, 1 - 5 mm, stubby handle
- 282354 DRANBS Straight Root Tip pick, standard handle
- 282386 DRANBS Tartar Removing Forceps
- 282387 DRANBS Extraction Forceps
- 282385 DRANBS Minnesota Retractor
- 282384 DRANBS Curette, subgingival
- 282383 DRANBS Jaquette Scaler, supragingival
- 282382 DRANBS Sickle Scaler
- 282380 DRANBS Explorer and Measuring Probe
- 282381 DRANBS Periosteal Elevator Double Ended

Art. No	Description
282358	DRANBS instrument set, 14 pieces, stubby handle



*Each instrument can be purchased individually*

**DRANBS Instrument Set, 5 Pieces, for Deciduous Teeth**

The set consist of:

- 282391 - 292393 DRANBS Curved Deciduous Elevator, S, M, L
- 282394 DRANBS Extraction Forceps, right angled
- 282390 DRANBS Periosteal Elevator, Molt, No 9

Art. No	Description
282359	DRANBS instrument set, 5 pieces, for deciduous teeth



*Each instrument can be purchased individually*

**DRANBS Feline Instrument Set (cats/small dogs), 8 Pieces**

The set consists of:

- 282347 - 282348 DRANBS Elevator, straight tip, 1 + 2 mm, stubby handle
- 282395 DRANBS Explorer and Measuring Probe, fine, cats/small dogs
- 282400 DRANBS Extraction Forceps, small, cats/small dogs
- 282396 DRANBS Periosteal Elevator, Double Ended, fine, cats/small dogs
- 282397 DRANBS Resorptive Lesion Probe
- 282398 DRANBS Sickle Scaler, fine, cats/small dogs
- 282399 DRANBS Curette, fine, cats/small dogs

Art. No	Description
282366	DRANBS feline instrument set (cats/small dogs), 8 pieces



*Each instrument can be purchased individually*

**DRANBS Elevator Set, Winged, 1-4 mm, Standard Handle**

The set consist of:

- 282330 DRANBS Elevator, winged, 1 mm, standard handle
- 282331 DRANBS Elevator, winged, 2 mm, standard handle
- 282332 DRANBS Elevator, winged, 3 mm, standard handle
- 282333 DRANBS Elevator, winged, 4 mm, standard handle

Art. No	Description
282360	DRANBS elevator set, winged, 1-4 mm, standard handle



*Each instrument can be purchased individually*

**DRANBS Elevator Set, Winged, 1-6 mm, Stubby Handle**

The set consist of:

- 282336 DRANBS Elevator, winged, 1 mm, stubby handle
- 282337 DRANBS Elevator, winged, 2 mm, stubby handle
- 282338 DRANBS Elevator, winged, 3 mm, stubby handle
- 282339 DRANBS Elevator, winged, 4 mm, stubby handle
- 282340 DRANBS Elevator, winged, 5 mm, stubby handle
- 282341 DRANBS Elevator, winged, 6 mm, stubby handle

Art. No	Description
282363	DRANBS elevator set, winged, 1-6 mm, stubby handle

*Each instrument can be purchased individually*



**DRANBS Luxator Set, 1-5 mm, Stubby Handle**

The set consist of:

- 282342 DRANBS Luxator, 1 mm, stubby handle
- 282343 DRANBS Luxator, 2 mm, stubby handle
- 282344 DRANBS Luxator, 3 mm, stubby handle
- 282345 DRANBS Luxator, 4 mm, stubby handle
- 282346 DRANBS Luxator, 5 mm, stubby handle

Art. No	Description
282364	DRANBS luxator set, 1-5 mm, stubby handle

*Each instrument can be purchased individually*



**DRANBS Serrated Luxator Set, 1 - 6 mm, Stubby Handle**

The set consist of:

- 282324 DRANBS Serrated Luxator 1 mm, stubby handle
- 282325 DRANBS Serrated Luxator 2 mm, stubby handle
- 282326 DRANBS Serrated Luxator 3 mm, stubby handle
- 282327 DRANBS Serrated Luxator 4 mm, stubby handle
- 282328 DRANBS Serrated Luxator 5 mm, stubby handle
- 282329 DRANBS Serrated Luxator 6 mm, stubby handle

Art. No	Description
282367	DRANBS serrated luxator set, 1-6 mm, stubby handle



*Each instrument can be purchased individually*

DRANBS Elevator Set, Winged, 1-6 mm, Stubby Handle, in Stainless Steel Case

The set contains:

- 282336 DRANBS Elevator, winged, 1 mm, stubby handle
- 282337 DRANBS Elevator, winged, 2 mm, stubby handle
- 282338 DRANBS Elevator, winged, 3 mm, stubby handle
- 282339 DRANBS Elevator, winged, 4 mm, stubby handle
- 282340 DRANBS Elevator, winged, 5 mm, stubby handle
- 282341 DRANBS Elevator, winged, 6 mm, stubby handle

Art. No	Description
282356	DRANBS elevator set, winged, 1-6 mm, stubby handle, in SS case
282369	DRANBS instrument case stainless steel, 16.5 x 14 x 3.7 cm



DRANBS Extraction Set, Winged Elevators with Stubby Handles, 12 pieces in Stainless Steel Case

The set consists of:

- 282336 – 282341 DRANBS Elevator winged, 1 – 6 mm, stubby handle
- 282381 DRANBS Periosteal Elevator double ended
- 282387 DRANBS Extraction Forceps
- 282408 DRANBS Needle Holder, Neivert
- 282349 DRANBS Tissue Forceps, Adson Brown
- 282409 DRANBS Scissor, Le Grange, serated
- 282350 DRANBS Scalpel blade handle, No 3

Art. No	Description
282365	DRANBS extraction set, winged elevators with stubby handles, 12 pieces in SS case

# Dranbs Impex

Mobile & Whatsapp  
+92 345 6815100

Email: [dental@dranbsimpex.com](mailto:dental@dranbsimpex.com)  
Web: [www.dranbsimpex.com](http://www.dranbsimpex.com)

Pakistan Production Office. G-237, Ugoki Model Town,  
Wazirabad Road, Sialkot. Pakistan